

Bat Dropping Sample Form



*Required Fields

Name*:

Company*:

Phone*:

Report Email*:

Invoice Email*:

☐

Same as report email

Invoice Address*:

PO/Reference:

New Customer/
Pay by Card:

☐

If you are a new customer **OR**
wish to pay by card please tick
this box.

Note: Analysis will NOT be started unless you have provided
us with a PO/Reference **or** ticked the Pay By Card option. For
samples in batches, each unique PO/Reference will get its
own report and invoice.

**For new customers, payment is required before results
are reported.**

Turnaround Time*:

☐

Standard

(10 Working Days)

☐

Fast Track

(5 Working Days)

See our website for current prices. Note: If no box is ticked, we
will assume a standard service required.

Sample 1)

Tube ID*:

Site Name:

O/S Reference:

Address:

Sampling Date:

Sample 2)

Tube ID*:

Site Name:

O/S Reference:

Address:

Sampling Date:

Sample 3)

Tube ID*:

Site Name:

O/S Reference:

Address:

Sampling Date:

Sample 4)

Tube ID*:

Site Name:

O/S Reference:

Address:

Sampling Date:

Sample 5)

Tube ID*:

Site Name:

O/S Reference:

Address:

Sampling Date:

Sample 6)

Tube ID*:

Site Name:

O/S Reference:

Address:

Sampling Date:

Instructions For Sample Collection

1. Conduct a thorough examination of the sampling location. Identify the dropping and reserve dropping required for analysis. Where possible, select droppings which are: **i) Out of direct sunlight ii) Complete and large iii) Fresh/recent.**
2. Put on the sterile gloves provided, using one of the spatulas provided, collect one whole dropping and scoop into one of the provided tubes. Collect another (or several) reserve droppings using the same spatula and place into the same tube. The reserve sample should be of a similar shape and size and likely be from the same species.
3. Close the tube lid ensuring it is screwed tightly and gently mix to ensure the droppings are placed in the preservative solution.
4. Write the sample I.D on the form provided correlating to the sample ID on the tube. (if tube is not pre-labelled ensure you label the tube using a marker pen or sticker).
5. Repeat steps above as required for up to another two species/sites using a separate tube/spatula for each sample.
6. Ensure you have filled in the form and included reserve samples where possible and then insert the completed form and samples in to the kit box and close the lid, sealing it with the included clear seals and tape (not provided). Post the freepost kit into a royal mail post box (postage is already paid).

For questions & support: **01332 292003** or **eDNA@SureScreen.com**

Have you seen our other DNA based services?

Mammals | Amphibians | Crayfish | Plants | Insects | Bird Sexing